

PRODUCT USE GUIDE FOR MEDISAN 70% HAND SANITIZING GEL D56

Product Description

A 70% alcohol hospital specification clear hand sanitizing gel.

Product Features

- Sanitizes the skin without the need for water.
- Contains emollients to prevent the drying of skin.
- Ideal for use in a medical environment, on entering food production areas, or before handling food.
- An alcohol based product for use in any situation where bactericidal control is important.
- Evaporates quickly and is non-tacky.
- Complies with BSEN 1500 Chemical Disinfectants and Antiseptics – Hygiene Handrub – Test Method.

Direction for Use

- Apply a small amount to the hands and rub in well, paying particular attention to finger nails and the area between the fingers.
- Allow hands to air dry.
- Important - the use of towels or tissue may lead to re-contamination with bacteria.
- Application to dry cracked or cut skin should be avoided.
- Thoroughly rinse the dispenser with warm water, pumping through the nozzle at least once a month to avoid blockage.

Storage

- Store in cool conditions away from direct sources of heat
- Avoid extremes of temperature
- Do not allow to freeze
- Ensure the cap is tightly fitted after use
- Keep out of reach of children

Shelf Life

- Shelf life 30 months, if stored as recommended.
- Refer to product label for further storage details

Environmental

- Do not allow large quantities of neat product to enter surface waterways.
- In the event of a large spillage, contact the local water authority or Environment Agency immediately.
- This product uses recyclable packaging. The detergents used in the formulation are biodegradable to EU guidelines under the detergents regulations (EC) No 648/2004.
- It contains no Phosphates.



Cleanol Group Ltd

Beaumont Road, Banbury, Oxon, OX16 1RB. UK.
Tel: +44 (0) 1295 251721 Fax: +44 (0) 1295 269561
Email: sales@cleanol.co.uk
www.cleanol.co.uk



Certificate No. 316



Certificate No. E3539

PRODUCT USE GUIDE FOR
MEDISAN 70% HAND SANITIZING GEL
D56

For further information please scan your QR Smart Reader here



Code Nos.

077123 2 x 5L
077137 6 x 500ml
077159 25 x 50ml

Data Sheet No. D56

Note: This product is only to be used for the purpose for which it was intended and in accordance with the specified instructions. This product should never be mixed with other chemical products and should be handled with care and stored carefully out of reach of children.



Cleenol Group Ltd

Beaumont Road, Banbury, Oxon, OX16 1RB. UK.
Tel: +44 (0) 1295 251721 Fax: +44 (0) 1295 269561
Email: sales@cleenol.co.uk
www.cleenol.co.uk



Certificate No. 316



Certificate No. E3539